

* Experiment Name :- SONOMETER

* Objective:- Sonometer with metallic wire, inbuilt AC supply, electromagnet, meter scale, slotted half kg weights, two knife edges.

* Theory and formula used:-

- Sonometer:- A sonometer is a device which consists of a thin wire stretched over two bridges that are usually mounted on a sound board and which is used to measure the vibration frequency, tension, density or diameter of the wire, or to verify relations between these quantities.
- Electromagnet:- An electromagnet is a type of magnet in which the magnetic field is produced by the flow of electric current. The magnetic field disappears when the current increases.
- Frequency of Ac mains:- The utility frequency or mains frequency is the frequency at which AC is transmitted from a power plant to the end user.

The natural frequency of the wire is given by,

$$f = \frac{1}{2l} \sqrt{\frac{T}{m}}$$

Teacher's Signature :

1905-01-02


20 min

20 min

2700°C (100) - sand turned

thin, very light glass turned - cut off
hot stem, temperate - glass is
soft, soft cut, give it flat bottom

- this shows if the melt



SONOMETER

is used for measuring tension of wire

where,

l = length of the sonometer wire between two knife edges when it is thrown into resonant vibrations.

T = Tension in the wire.

m = mass per unit length of the wire.

The frequency of the AC mains is given by,

$$n = \frac{f}{2} = \frac{1}{4l} \sqrt{\frac{T}{m}}$$

* Observation:- Mass per unit length of the wire (m) = 0.0056 kg/m.

S.No	load in (kg)	Tension in Newton $T = (w \times 9.8)$	Resonating Position of bridges increasing (l_1)	length of the wire(cm)	Frequency
			Position of bridges decreasing (l_2)	Mean resonating length (l)	$f = \frac{1}{2l} \sqrt{\frac{T}{m}}$
	2.1 kg	20.58	28.1	26.5	27.5
	1.6 kg	15.68	27.5	25.5	26.5
	1.1 kg	10.78	27.1	27.45	27.05
	0.6 kg	5.88	13.2	13.66	12.6
					123.1

Teacher's Signature : _____

* Calculation :-

$$\textcircled{1} \quad f = \frac{1}{2l} \sqrt{\frac{T}{m}}$$

$$f = \frac{1}{2 \times 27.5} \sqrt{\frac{20.50}{0.0056}}$$

$$f \Rightarrow 112.25 \text{ Hz}$$

$$\textcircled{2} \quad f = \frac{1}{2l} \sqrt{\frac{T}{m}}$$

$$f = \frac{1}{2 \times 26.5} \sqrt{\frac{15.68}{0.0056}}$$

$$\Rightarrow 99.83 \text{ Hz}$$

$$\textcircled{3} \quad f = \frac{1}{2l} \sqrt{\frac{T}{m}}$$

$$f \Rightarrow \frac{1}{2 \times 27.05} \sqrt{\frac{10.78}{0.0056}}$$

$$f = 81.16 \text{ Hz}$$

$$\textcircled{4} \quad f = \frac{1}{2l} \sqrt{\frac{T}{m}}$$

$$f = \frac{1}{2 \times 12.6} \sqrt{\frac{5.88}{0.0056}}$$

$$f = 123.12 \text{ Hz}$$

Teacher's Signature : _____

$$\text{Mean} = \frac{112.25 + 99.83 + 101.16 + 123.12}{4}$$

$$= 104.152 \text{ Hz}$$

frequency of AC mains $\Rightarrow n = \frac{f}{2}$

$$\Rightarrow \frac{104.152}{2}$$

$$\Rightarrow 52.076$$

* Result : The frequency of AC mains = 52.076 Hz
 Standard Result \rightarrow The frequency of AC Mains = 50 Hz

* Percentage error :-

$$\Rightarrow \left| \frac{\text{standard value} - \text{calculated value}}{\text{standard value}} \right| \times 100$$


$$\Rightarrow \left| \frac{50 - 52.06}{50} \right| \times 100$$

$$\Rightarrow 0.04 \times 100$$

$$\Rightarrow 4\%$$

* Precaution and source of error :-

i) The wire should be of a uniform area of cross-section free from kinks and should be light.

 Teacher's Signature : _____


- ii. The observation should start with maximum distance between the two knife edges.
- iii. The weight of hanger should always be included in the load.
- iv. The pulley should be free from friction.



* SOLARCELL :-

- Objective :- To draw the characteristics of solar cell and to estimate fill factor (FF) of solar cell.
- Apparatus Required :- Solar cell, rheostat, ammeter, voltmeter, illumination source, bulb and connecting wires.
- Circuit Diagram :-
⇒ next page.
- Theory and formula used :- solar cell is basically a two terminal p-n junction device designed to absorb photon absorption through the electrical signal or power in the external circuits. Therefore it is necessary to discuss the physics of semiconductor p-n junction diode, which converts the optical energy into electrical signals.

Photo voltaic system converts sunlight directly into electrical energy. The backbone of this technology is semi-conducting materials such as silicon. A typical solar cell consists of two



* Set up for determining characteristics of a solar cell

differently doped semiconductor p-n junction diode, which converts the optical energy into the host material. Starting out with a pure semiconductor crystal. This is achieved by substituting some of the atoms in the crystal lattice with elements that have one more or less valence electrons than the host material.

~~semiconductor~~ elements have four valence electrons all of which are used for bonding in the crystal lattice. If doping is controlled introduction of impurities into the host material. If doping material has five electrons all of which are used for bonding in the crystal lattice. These 'free' atoms can move about easily in lattice and are responsible for an increase in conductivity.

Since they have a negative charge the doping material doped in this way is called n-type semiconductor. Since in case of free charge carriers are positive this kind of semiconductor is said to be of p-type. When a p-type semiconductor is joined to an n-type semiconductor, a p-n junction is created.

* Characteristics of a solar cell :-


The usual voltage from solar cells depends on the semiconductor material. In silicon it amounts to approximately 0.5 V. Terminal.

voltage is only weakly dependent on light radiation, while the current intensity increases with higher luminosity.

* Fill factor (FF) :-

FF percentage measures the "squareness" of the I-V curve. It states the degree to which the voltage at the maximum power point (V_{mp}) matches the open-circuit voltage (V_{oc}) and that the current at the maximum power point (I_{mp}). This relation is given by:-

$$FF = \frac{V_m I_m}{V_{sc} I_{sc}}$$




* Efficiency of solar cell :-

The efficiency of the solar cell is the ratio of produced electrical power (P_{out}) and the incident radiant power (P_{in}).

Efficiency of solar cell,

$$\eta = \frac{P_{out}}{P_{in}}$$




Observation Table :-

SNo.	Voltage, V (Volts)	Current, I (mA)	Power calculated $P = VI$ (watt)	$P = V \cdot I$ (watt) $P = VI / 100$
1.	1	119.6	119.6	1.196
2-	2	117.6	235.2	2.352
3-	3	114.8	344.4	3.444
4-	4	112.1	448.4	4.484
5-	5	110.2	551.0	5.510
6-	6	109.1	684.6	6.846
7.	7	109.4	779.8	7.798
8-	8	108.9	871.2	$\Rightarrow 8.712$
9-	9	94.8	853.3	8.533
10-	9.90	14.1	139.59	1.395
12-	0.12	124.6	14.592	0.145

*** Procedure:-**

1. Take the solar Energy trainer along with solar panel.
2. Place the solar panel in the stand and adjust the panel at an angle of about 45° with the ground. Direct the sunlight straight at the solar panel (angle of 90°)
3. With the DB 15 - connector connect the solar energy trainer with solar panel. Then wait for 1 minute to avoid errors due to temperature fluctuations.
4. Set the potentiometer to maximum resistance, i.e., at fully clockwise position and measure and record its resistance into the observation table.
5. Connect solar cell.
6. Connect positive terminal of solar cell to P₁ terminal of potentiometer.
7. Connect other end of potentiometer i.e., P₂ to positive terminal of ammeter.

Teacher's Signature : _____



at 29°C extrapolated for base resistance = 0.01

8E Connect negative terminal of ammeter to negative terminal of solar cell.

9E Now connect the positive terminal to P₁ and negative terminal to P₂.

10E Record the values of corresponding voltage and current into observation table.

11E. The Fill Factor (FF) is defined by the ratio

$$FF = \frac{V_m I_m}{V_{sc} I_{sc}}$$

$$\Rightarrow \frac{0 \times 108.9}{10.49 \times 121}$$

$$\Rightarrow \frac{0.712}{1269.29}$$

$$\Rightarrow 0.6863$$

* Calculation:-

1E. Plot the IV characteristic of solar cell.


2E. Plot the power (VI) vs V and determine the maximum power.

Ans 0.112121

Teacher's Signature : _____

3- Determining the Fill factor(FF).


$$FF = \frac{V_m I_m}{V_{sc} I_{sc}}$$



* Results:-

1) The I-V characteristic was drawn for given solar cell.

2) Fill factor(FF) = 0.6863




* Precautions:-


1) Make sure that all the connections are tight.

2) Wait for few minutes to avoid errors due to temperature fluctuations.


3) To measure the resistance of potentiometer at any position.



Teacher's Signature : _____



(a) Forward-bias



(b) Reverse-bias

P-n Junction:

- Objective:- To draw the (V-I) characteristics of P-n junction diode and to estimate the dynamic resistance.
- Apparatus required:- Diode trainer kit, Power supply voltmeter, ammeter and connecting wires.
- Theory used:- An equation describes the exact current through a diode, given the voltage dropped across the junction the temperature of the junction and several physical constants. It is commonly known as the diode equation.

$$I_D = I_s \{ \exp (eV_D / k_B T) - 1 \}$$

where I_D is diode current in amperes and I_s is total saturation current in amperes.

V_D = Voltage applied across diode in volts.
It is positive for forward bias and negative for reverse bias.


k_B = Boltzmann's constant ($1.38 \times 10^{-23} \text{ J/Kelvin}$)

T = Junction temperature in Kelvin.

Such that, $\exp (eV_D / k_B T) \gg 1$ then

$$I_D = I_s \exp (eV_D / k_B T)$$

Teacher's Signature :



This shows that in forward bias, current increase exponentially as shown.
(V-I) graph in figure.

Observation :-

S.N.O.	Forward bias (1)		Reverse bias (2)	
	Diode voltage V_D	Diode current (mA)	Diode voltage •	Diode current (mA)
1.	0.1	0	1	2.03
2.	0.2	0	2	3.49
3.	0.3	0	3	4.83
4.	0.4	0.1	4	6.10
5.	0.5	2.1	5	7.49
6.	0.6	30.4	6	8.83
7.	0.7	58.6	7	10.28
8.	0.8	128.4	8	11.68
9.	0.9	181.2	9	13.18

calculation :- Plot a graph between V and I.


Find the static and dynamic resistance.

Static resistance :-

$$R_D = V_D / I_D$$

$$\Rightarrow \frac{0.6}{30.4} \Rightarrow 0.019 \Omega$$


Teacher's Signature : _____



- Result:- The I-V characteristics of the diode are shown in the forward and reverse bias. The static and dynamic resistance are :- 0.01952

- Precautions:-

- ① Make sure that all the connections are tight.
- ② The voltage should not exceed the specified breakdown voltage in reverse bias.



Hall Effect

- Objective:- To study the hall effect and to determine the hall co-efficient, charge carrier concentration and carrier mobility of given semiconductor material using hall set up.
- Apparatus Required:- Gauss and Tesla meter, measurement unit, constant current power supply, electromagnet, Hall probe, In-As probe, etc.
- Theory and formula used:- If an electric current flow through a conductor in a magnetic field, the magnetic field, the magnetic exerts transverse force on the moving charge carriers which tends to push them to one side of the conductor. A build up of charge at the sides of the conductor will balance this magnetic influence, producing a measurable voltage between the two sides of the conductor. The presence of this measurement transverse voltage is called the Hall effect after 1879.

The Hall coefficient is given by,

$$R_H = \frac{V_H d}{IB} \quad \text{--- (1)}$$

Teacher's Signature : _____

Units, $V_H \rightarrow$ Hall voltage $I \rightarrow$ Probe current $d \rightarrow$ Thickness of sample $B \rightarrow$ applied magnetic field,

- Observations :-

S.NO.	I_B (A)	Magnetic field (B^+) from one side of probe)	Magnetic (B^-) field [from another side of probe]	mean \Rightarrow $(B^+ - B^-)/2$ B (Tesla)
1.	1.0	1.62	-1.62	1.62
2.	1.5	1.91	-1.91	1.91
3.	2.0	2.36	-2.36	2.36
4.	2.5	2.63	-2.63	2.63
5.	3.0	3.02	-3.20	3.11
6.	3.5	3.01	-3.21	3.11
7.	4.0	3.24	-3.40	3.32
8.	4.5	3.40	-3.56	3.48

Teacher's Signature : _____

Expt. No. _____

Date _____

Page No. _____

S.No.	Probe current (mA)	Zero field potential at set voltage V_0 , (mv)	Observed Voltage for one side (V_{t_n})	Observed voltage for another side (V_{H_n})	Actual Hall Voltage $(V_H = V_H^+ - V_0)$	Actual Hall Voltage $(V = V_H^- - V_0)$	Measured Voltage, mean voltage $V_H - V$
1.	1.0	17.6	30.0	5.3	12.4	-12.3	12.35
2.	1.5	26.9	43.8	8.7	16.9	-18.2	17.55
3.	2.0	36.4	56.4	14.1	20.0	-22.3	21.15
4.	2.5	43.1	70.3	14.7	27.2	-28.4	20.95
5.	3.0	50.1	82.4	19.0	32.3	-31.1	31.70
6.	3.5	58.1	94.3	10.1	36.2	-40.0	38.10
7.	4.0	65.	104.4	22.6	31.4	-42.6	37.00
8.	4.5	73.0	122.5	24.4	49.5	-48.6	49.05
9.	5.0	78.5	130.0	26.0	51.5	-52.5	52.00

• Formulas :-

$$R_H = \frac{V_H * t}{I * B}$$

where, $V_H \Rightarrow$ Hall voltage

B = applied magnetic.

$t \Rightarrow$ Thickness

R_H = Hall co-efficient

$I \Rightarrow$ current

$$\boxed{R_H = \frac{1}{n e}}$$

where, $e \Rightarrow$ electronic charge

$n \Rightarrow$ number density of electrons in conductor.

Teacher's Signature : _____

The charge 'e' carrier density 'n' can be calculated,

$$n = \frac{1}{R_H e}$$

$$1/R_H e$$

We know that carrier mobility is the velocity of current carrying particles per unit electric field and it is related with Hall coefficient as -

$$m = \sigma R_H$$

$$\text{where, } \sigma = \frac{1}{\rho} \text{ or conductivity}$$

$\rho \rightarrow$ resistivity.

$$m = \frac{1}{\rho} R_H$$

$$\Rightarrow \frac{1}{0.47} \times 1000 \text{ esu}$$


$$\Rightarrow 3.33 \times 10^{-4} \text{ cm}^2/\text{volt sec.}$$

- Results \rightarrow The Hall co-efficient (R_H) = $1.569 \times 10^{-1} \text{ cm}^3/\text{coulomb}$.
 The value of charge carrier density = $3.98 \times 10^{19} \text{ cm}^{-3}$
 The value of charge mobility (m) = $3.33 \times 10^{-4} \text{ cm}^2/\text{volt sec.}$

15.12.21

Teacher's Signature : _____

Winkel θ ist der Winkel zwischen dem Vektor zur Richtung des einfallenden Lichtstrahls und dem Vektor zur Richtung des reflektierten Lichtstrahls.



18.08.2018

Experiment →

Diffraction Grating

- Objective:- To find the wavelength of monochromatic light with the help of a plane transmission diffraction grating and spectrometer.
- Apparatus:- A spectrometer, a spirit level, a ~~sq~~ source of monochromatic light (sodium lamp), diffraction grating with clamping arrangement.
- Formula used:- The wavelength ' λ ' of any spectral line can be calculated by the formula.

$$(a+b) \sin\theta = n\lambda$$

$$\lambda = \frac{(a+b) \sin\theta}{n}$$

where $a+b$ = grating element

θ = Angle of diffraction

n = order of spectrum

Theory :- A plan diffraction grating consists of an optically plan glass plate on which large no. of equidistant parallel lines are ruled. These lines divide the glass plate into opacities and transparencies, the thickness of which are of the order of wavelength of the visible light. The region where a line is drawn becomes opaque and the space between the two lines is transparent.

When a parallel beam of monochromatic light is incident normally on the grating. If it suffers diffraction, the transmitted light gives rise to primary maxima in certain directions given by the relation,

$$(a+b) \sin\theta = n\lambda$$

where $a \rightarrow$ width of transparency
 $b \rightarrow$ is that of opacity. θ is the angle of diffraction for n^{th} order maxima and λ is the wave length of light.

For the first order spectrum, $n=1$

$$\text{Therefore, } (a+b) \sin\theta_1 = \lambda$$

for the second order spectrum, $n=2$

Teacher's Signature : _____

Order of vernier spectrum scale	No.	Reading of Telescope from left side Spectrum			Reading of Telescope from right side Spectrum			Angle $\theta = \frac{(a-b)}{4}$	Mean Angle θ	Angle θ
		M.S. reading	V.S. reading	Total reading	M.S. reading	V.S. reading	Total reading			
first	V_1	59.5	9	59.65	94	10	94.16	34.51	34.4	17.24
	V_2	240	4	240.06	274.5	1	274.5	34.44		
Second	V_1	64	5	64.003	89	20	89.33	25.27	25.16	12.58
	V_2	244.5	0	244.5	269.5	3	269.5	25.05		

and $(a+b) \sin\theta_2 = 27$

The number of lines N per inch are marked on the grating element $(a+b)$ is given by:-

$$(a+b) = 2.54 \text{ cm.}$$

Spectrometer :- An spectrometer is basically an instrument for measuring the angular deviation of light ray. A prism produces this angular deviation of light ray or grating, which is wavelength dependent. The emergent light from these component is dispersed into a spectrum in which wavelength is a function of angle. Spectrometer consists of mainly three parts : collimator, Prism table and ~~left telescope~~.

Observations:-

(A) - Determination of least count of spectrometer -
Reading of one smallest division of M.S. ' a ' = 0.5 degree
Total no. of division on vernier scale ' b ' = 30

$$\text{Least count of spectrometer} = \frac{a}{b} \rightarrow \frac{0.5}{30} \text{ degree}$$

\Rightarrow 1 minute

Teacher's Signature : _____

(B) Determination of grating element of diffraction grating →

No. of lines per inch on the grating $N = 300$

$$\therefore \text{Grating element } (a+b) = \frac{2.54}{N} \text{ cm}$$

$$= \underline{1.693 \times 10^{-4} \text{ cm}}$$

(C) Determination of angle of diffraction →

Calculation :-

First order spectrum:-

wavelength of spectral lines of first order $\star (n=1)$

$$\lambda_1 = (a+b) \sin \theta$$

$$\lambda_1 = \underline{4662 \text{ Å}}$$

Second order spectrum:-

wavelength of spectral lines of second order $(n=2)$

$$\lambda_{12} = \frac{(a+b) \sin \theta_2}{2}$$

$$\lambda_2 = \underline{7036 \text{ Å}}$$

$$\text{Mean of } \lambda = \frac{(\lambda_1 + \lambda_2)}{2} : \Rightarrow 5852 \text{ Å}$$

Teacher's Signature : _____

Result :-

The wavelength of given monochromatic light
 $\lambda \rightarrow 3673 \text{ \AA}$


Percentage error:-

$$\% \text{ Error} = \frac{(\text{Standard value} - \text{Experimental value})}{\text{Standard value}} \times 100$$

Precautions:-

The slit should be narrow

The diffraction grating should be parallel to the slit.



Experiment - 6

He-Ne Laser


Objective :- To determine the wavelength of He-Ne laser source using diffraction grating.

Apparatus :- HeNe laser source, Diffraction grating with its mount, Optical bench, a screen with millimeter scale.


Theory :- Diffraction is the bending of wave around an obstacle or small aperture whose size is comparable to the wavelength of light if source and screen are far away from each other than Fresnel diffraction occurs.

A plane diffraction grating consists of an optically plane glass plate on which a number of equidistant lines are ruled. These lines divide the glass plate into opaque and transparent places, the width of which is of the order of wavelength of visible light. The region where a line is drawn become opaque whereas the space between two lines is transparent and acts as

Teacher's Signature : _____



Experimental setup for
diffraction



slit. The diffraction maxima satisfying the grating condition.

$$(a+b) \sin \theta_n = n\lambda$$

where $(a+b)$ \rightarrow grating element

a' \rightarrow separation between two sets.

'b' \rightarrow slit width.

A. $d = (a+b)$

$$d \sin \theta_n = n\lambda$$

Working formula:- If a parallel beam of monochromatic light is incident normally on the plane diffraction grating, maxima & minima occurs due to diffraction. The condition for diffraction maxima is \rightarrow

$$d \sin \theta_n = n\lambda \quad \text{--- (1)}$$

where d is known as grating element,

θ_n is the angle of diffraction of n^{th} maximum.

n is order of spectrum.

λ is the wavelength of incident light, again

$$\tan \theta = \frac{y}{L}$$

Hence $\sin \theta = \frac{y}{\sqrt{y^2 + L^2}} \quad \text{--- (2)}$

Teacher's Signature : _____

* Observations →

No. of lines per centimeter on grating, $N = 3000/\text{cm}$

$$\text{Grating element, } d = \frac{1}{N} = \frac{1}{3000} \text{ cm}$$

S. No.	Distance b/w screen and grating L (cm)	order of maxima (n)	Distance of maxima from centre			$\lambda = \frac{d \cdot y}{\sqrt{y^2 + L^2}}$ (in Å)	λ_{mean} (Å)
			Left side (Y ₁)	Right side (Y ₂)	Mean, $y = \frac{(Y_1 + Y_2)}{2}$		
1.	10	1	2	2	$\frac{2+2}{2} = 2$	6542.36	
2.	15	1	3	3	$\frac{3+3}{2} = 3$	6540.22	
3.	20	1	3.8	4.6	$\frac{3.8+4.6}{2} = 4.2$	6381.93	6441.14 Å
4.	25	1	4.9	4.9	$\frac{4.9+4.9}{2} = 4.9$	6412.77	
5.	30	1	5.8	5.8	$\frac{5.8+5.8}{2} = 5.8$	6328.42	

using equation(2)

$$d = \frac{y}{\sqrt{y^2 + L^2}} = n\lambda$$

Hence, $\lambda = \frac{dy}{n\sqrt{y^2 + L^2}}$

$$\lambda = \frac{1}{N} ; N \text{ is number of slits/cm.}$$

Calculations :-

$$\lambda = \frac{dy}{n\sqrt{y^2 + L^2}}$$

$$\lambda_{\text{mean}} = \frac{\lambda_1 + \lambda_2 + \lambda_3 + \dots + \lambda_n}{n}$$

$$\Rightarrow 6.4414 \times 10^{-5} \text{ cm}$$

$$\Rightarrow 6441.14 \text{ Å}$$

Result :-

The wavelength of given laser source,

$$\lambda = 6441.14 \text{ Å}$$

The standard value of Hene laser source
is 6328 Å

Percentage error :-


$$\text{percentage error} = \frac{6441.14 - 6328}{6328} \times 100$$

$$\Rightarrow \frac{113.14}{6328} \times 100$$

$$\Rightarrow 1.78\%$$

Teacher's Signature : _____

Circuit Diagram



Experiment -07Planck's Constant

- Objective :- To determine planck's constant using light emitting diode.
- Apparatus used :- Trainer kit, super bright LEDs and variable power supply (0-5) Volts.
- Objective :-
Theory & formula used →
In this experiment, we are using LED to measure planck's constant. LED's are semiconductor diode that emits electromagnetic radiation in optical and near optical frequencies, when operated in forward bias, above a minimum threshold voltage. In this condition, an electron hole pair is created in diode & thus current starts flowing. Above the threshold value, the current increases exponentially with voltage. A quantum of energy is required to create an electron hole pair and this energy is released when an electron & a hole recombine.

Teacher's Signature : _____

* Observation Table →

S. No.	Colour	Wavelength (A)	$\frac{1}{\lambda}$ ($m^{-1} \times 10^{-6}$)	V _o (Volts)
1 -	Blue	470 nm	2.127	2.23
2 -	Green	510 nm	1.961	2.20
3 -	Yellow	570 nm	1.754	1.72
4 -	Orange	600 nm	1.666	1.56
5 -	Red.	700 nm	1.428	1.43

The energy of photon is related to its wavelength. a:-

$$E = \frac{hc}{\lambda} = h\nu$$

$$\frac{hc}{\lambda} = ev_0 \quad \text{or} \quad h = \frac{ev_0\lambda}{c}$$

- Results :- Planck's constant is $7.066 \times 10^{-34} \text{ Js}$
Standard value is $6.625 \times 10^{-34} \text{ Jc}$

$$\% \text{ error} = \left| \frac{(6.625 - 7.066) \times 10^{-34}}{6.625 \times 10^{-34}} \times 100 \right|$$

$$\Rightarrow \left| -0.441 \times 100 \right|$$

$$\rightarrow 6.65\%$$

- Sources of error →
- Make sure that all connections are tight
- Note the threshold voltage carefully, when the LEDs just start flow.
- Handle the equipment carefully.
- Make sure that LED is in the forward bias